

# INNOVATION THAT MAKES SENSE

## Integrated Flowmeter...

Gas Specific Labeling and Fittings  
Medical Oxygen and Medical Air



Supply Outlet  
50 PSI

Adjustable Knob for flow control  
Ranges available:  
0-8 LPM, 0-15 LPM, 0-25 LPM

Flow Indicator

Flowmeter Port

- ✓ 2 outlets in 1
- ✓ Easy to service
- ✓ Retrofit existing outlets
- ✓ No more lost, chained, or broken Flowmeters
- ✓ Safest, most convenient way to dispense Medical Oxygen & Air

Amico Integrated Flowmeter outlet combines a Dial Flowmeter and a Medical Gas Outlet into a single compact design. The right side port connection has a control knob which allows the user to adjust the flow setting, while the left side port is a direct connection to the gas supply.

The Amico Integrated Flowmeter saves space, reduces construction costs, and ensures healthcare providers always have a Flowmeter when they need one.

Part Configurator

## S-FI-CXXYZ-WQ-S

**Flow range:**

8 LPM: 08  
 15 LPM: 15  
 25 LPM: 25

**Colour Code:**

USA: U  
 ISO: I

**Gas Type:**

Oxygen: O  
 Medical Air: A

**Line Fitting:**

DISS Male  
 with Check valve: D  
 Ohmeda Coupler: O  
 Chemetron Coupler: C

**Flowmeter Fitting:**

DISS Male  
 with Check valve: D

**Amico Retro-Fit Conversion Compatibility:**

Diamond 1 & Diamond 2/3 (Complete Conversion):	A
Diamond 2/3 (Plug and Play):	B
Puritan Bennett Series 73, Series A & Litharge NCG/Chemetron Series 500, 560 Oxequip/Medstar Series 27 (Smooth Bore):	C
NCG/Chemetron Series 148, 400, 460 Oxequip/Medstar Series 7 :	Leave Blank
(Regular) Amico Rough-in:	Leave Blank

NOTE: Rough-in conversions, Filler Frames and Trim plates are not included in these part numbers

represented by: