

MEDICAL GAS PIPELINE

RETRO-FIT

PRODUCTS

DON'T CUT
YOUR PIPES!

DON'T BREAK
YOUR WALLS!

INCREASE FLOW

LEAKAGE
PROBLEMS?

CAN'T GET PARTS?

AGING PARTS

DOING
RENOVATIONS?

UPGRADE
TECHNOLOGY

SAVE ON REPAIRS




ALARM RETRO-FIT



Retro-Man says: “Medical Gas Alarms are designed to update and convert your old alarms to a new Amico Area, Master or Combination Alarms.”

Many companies have discontinued their older Alarms, therefore spare parts to fix these aging systems are a rare commodity. Amico’s “Retro-Fit” is the solution you’ve been looking for. By retro-fitting a facility you can replace the existing Alarm without cutting the pipeline or breaking into the existing wall. Amico’s Retro-Fits are designed to use the existing back box of the previous product. This allows you to maintain the original configuration, while improving both the functionality and aesthetic appearance of the Alarm.



Chemetron Digital Alarm



Amico Alarm



Chemetron Analog Alarm




Amico Alarm



Ohio Alarm




Amico Alarm



Ohmeda Alarm



Amico Alarm



Medæx Integrated Alarm



Amico Alarm



Oxequip Alarm




Amico Alarm



Puritan Bennett Analog Alarm



Amico Alarm



Puritan Bennett Digital Alarm



Amico Alarm

OUTLET RETRO-FIT



Retro-Man says: “Our Medical Gas Outlets are designed to update and convert your old Outlets.”

The Retro-Fit process allows you to update and convert your old Medical Gas Outlets to new Amico Outlets. This can be done without cutting the pipeline or breaking into the existing wall. Amico’s Retro-Fits are designed to use the existing back body of the previous product. This allows you to maintain the original configuration, while improving both the functionality and aesthetic appearance of the Outlet.



Diamond 1 Outlet



Amico Ohmeda Compatible or DISS Outlet



Diamond 2 Outlet



Amico Ohmeda Compatible or DISS Outlet



Diamond 3 Outlet



Amico DISS or Ohmeda Compatible Outlet



NCG Outlet




Amico Ohmeda, Chemetron, Puritan Bennett Compatible or DISS Outlet



Oxequip Outlet



Amico Ohmeda, Chemetron, Puritan Bennett Compatible or DISS Outlet



Puritan Bennett Outlets



Amico Ohmeda, Puritan Bennett Compatible or DISS Outlet

GAS CONTROL PANEL RETRO-FIT



Retro-Man says: “Our Medical Gas Control Panels are designed to update and convert your old Gas Control Panels.”

The Retro-Fit process allows you to update and convert your old Medical Gas Control Panel to a new Amico Gas Control Panel. This can be done without cutting the pipeline or breaking into the existing wall. Amico’s Retro-Fits are designed to use the existing back box of the previous product. This allows you to maintain the original configuration, while improving both the functionality and aesthetic appearance of your Gas Control Panel.



Puritan Bennett II Gas Control Panel



Amico Gas Control Panel



Chemetron Gas Control Panel



Amico Gas Control Panel



Ohmeda Gas Control Panel



Amico Gas Control Panel



Oxequip Gas Control Panel



Amico Gas Control Panel


THE ULTIMATE RETRO-FIT



Retro-Man says: “The Ultimate Retro-Fit will allow you to update and/or add new products such as: Alarms, Outlets, Gas Control Panels and Valve Boxes to existing locations. This design will surface mount to any wall, like a picture. Installation is quick and easy. Just pipe in from above the ceiling or behind the wall.”



DISTRIBUTED BY:



85 Fulton Way, Richmond Hill, ON L4B 2N4 Canada • Tel: 1-877-462-6426 • Fax: 905-764-0862
31 Sprague Avenue, Amityville, NY 11701, USA • sales@amico.com • www.amico.com