





AHM monitor arm



Features:

- Parallel linkage provides consistent viewing angle at any position
- 17" (43.2 cm) of height adjustment
- 180° degrees side to side motion
- Internal gas spring provides counter balance
- Lock arm in any position
- 30 lbs (14kg) arm rating — heavy duty option available for loads up to 40 lbs (18.2 kg)
- Can be mounted directly to any equipment rail or wall channel
- Integrated cable management
- Color options available (white is standard)



Monitor Compatibility

Patient Monitors	LCD Monitors
• Draeger	• Toshiba
• Spacelabs	• Sony
• GE	• Samsung
• Nihon Kohden	• Gateway
• Datex Ohmeda	• Advantech
• Welch Allyn	• DELL
• Datascope	• NEC
• Philips	• LG
Other Manufacturers Options Available	

Please specify model when ordering.

Channel Mount Options

Channel to **Wall**



- Channel mounts directly to the wall

Channel to **Equipment Rail**



- Channel mounts directly to equipment rail
- Channel can be mounted anywhere along length of rail

Channel to **Power Column**



- AHM mounts directly to existing monitor channel on power column or flatwall

Channel to **Headwall**



- Channel mounts directly to existing rail on headwall
- Channel can be mounted anywhere along length of headwall

Direct Rail Mount Options

Benefits of **Direct Rail Mounting**

- Eliminates the need for a monitor channel
- Less clinical look
- Increases flexibility for location of monitor
- 30 lbs (14 kg) maximum weight rating
- Compatible with various headwall, flatwall, equipment booms and rail systems

Direct to **Equipment Rail**



- Mounts directly to equipment rail



Direct to **Power Column**



- Mounts directly to existing vertical rail on power column or flatwall



Direct to **Headwall**



- Mounts directly to existing horizontal rail on headwall



Additional **AHM Models**



LCD AHM Monitor Arm

- 18" (45.7cm) extension
- 30 lbs (13.60kg) arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing

Part Number: AHM-LCD-05



AHM Monitor Arm

- 18" (45.7cm) extension
- 30 lbs (13.60kg) arm rating
- $\pm 30^\circ$ degrees tilt

Part Number: AHM-050



LCD AHM Monitor Arm with 9" extension

- 27" (68.6cm) extension
- 30 lbs (13.60kg) arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing

Part Number: AHM-LCD-35



AHM Monitor Arm with 9" extension

- 27" (68.6cm) extension
- 30 lbs (13.60kg) arm rating
- $\pm 30^\circ$ degrees tilt

Part Number: AHM-350



LCD AHM Monitor Arm with 12" extension

- 30" (76.2cm) extension
- 30 lbs (13.60kg) arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing


Part Number: AHM-LCD-15



AHM Monitor Arm with 12" extension

- 30" (76.2cm) extension
- 30 lbs (13.60kg) arm rating
- $\pm 30^\circ$ degrees tilt

Part Number: AHM-150



LCD AHM Monitor Arm with 16" extension

- 34" (86.4cm) extension
- 30 lbs (13.60kg) arm rating
- VESA 75/100 plates included
- Portrait and landscape viewing

Part Number: AHM-LCD-55



AHM Monitor Arm with 16" extension

- 34" (86.4cm) extension
- 30 lbs (13.60kg) arm rating
- $\pm 30^\circ$ degrees tilt

Part Number: AHM-550

AHM Accessories

Power Bar



- Hospital grade, 4 power outlets
- 120 VAC, 60Hz

Part Number: POWER4-MON

Power Holder



- Holds most power adapters for PCs, patient monitors and LCDs

Part Number: PWRHLD-MON

CPU Mount



- Conveniently mounts CPU to wall channel ranging in sizes from 1.5" (3.8cm) to 9.5" (45.7cm)
- Call with dimension of CPU width

Call for part number.

Cord Hooks



- Manage cables, cords and other hanging wires
- Attach to the bottom of any monitor arm

Part Number: CRDHLD

Monitor Channel



- Standard wall channel 19" (48.3cm)
- Seismic options are available and custom lengths are welcome

Part Number: MRS-019


Channel Covers



- Easy way to hide wires
- Available in patterns and colors to match any decor

Call for part number.

Small Storage Basket



- 10"W x 10"D x 4" H white wire basket
- Can be sterilized

Part Number: 101004-MON-BB


BP Cuff Basket



- Conveniently holds single BP Cuffs to wall channel
- Can be sterilized

Part Number: CUFBASK-MON


Medium Storage Basket



- 10"W x 06"D x 09" H white wire basket
- Can be sterilized

Part Number: 100609-MON-BB

Large Storage Basket




- 13"W x 09"D x 07" H white wire basket
- Can be sterilized

Part Number: 130907-MON-BB

Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more mounting solutions.



Distributed By:



SP-AHM063