


Counter-Top mounts



Part Number: CTM-12-A00

Features:

- 8", 12", 18" and 24" height profiles and custom lengths available
- Internal raceway for cable management
- Double sided profile allows for any application
- Holds up to 100lbs (45kg) flush to column and 50lbs (23kg) 24" away
- 8" x 8" base plate, plus an additional backing plate for secure installation
- Two (2) side raceway covers and two channel covers (field cut to desired lengths)
- Clear anodized aluminum finish with optional colors

Dual Channel Counter-Top Mounting Applications

AHM Mounted



AHM/Extension Arm Mounted



SSM Mounted



Height Profiles

- 8" Part Number: CTM-08-A00
- 12" Part Number: CTM-12-A00
- 18" Part Number: CTM-18-A00
- 24" Part Number: CTM-24-A00

Optional Colors



Table-Top mounts



Part Number: TTM-12-450-A00

Features:

















- 3", 5", 9", 12", and 16" height profiles
- Compatible with flat panel displays and patient monitors
- Base plate, plus an additional backing plate for secure installation
- 1 1/4" Diameter clear anodized aluminum pole
- Adjustable tilt head allows multiple viewing angles
- Holds up to 30lbs
- Available in custom colors and different lengths

Height Profiles

3" Part Number: TTM-03-450-A00
 5" Part Number: TTM-05-450-A00
 9" Part Number: TTM-09-450-A00
 12" Part Number: TTM-12-450-A00
 16" Part Number: TTM-16-450-A00

Table -Top Mount Part Matrix

TTM-WW-45Y-ZZZ

 03	 09	 16	 0	 6	Standard Clear Anodized  A00 Pewter  A01 Forest Green  A02 Blue  A03 Red Maroon  A04
 05	 12		 2	 7	
			 3	 8	

Distributed By:

Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more mounting solutions.

