

Equipment roll stand



Features:

- Load capacity up to 25 lbs (5.5kg)
- Overall height 42" (106.7cm)
- 3" (7.6cm) single wheel casters (2 locking, 3 non-locking)
- 1 or 2 foot horizontal rail options available
- Custom color options available (grey is standard)
- MRI safe option available
- Sturdy 25" aluminum base
- Variety of accessories available

ARS Equipment Roll Stand



- Amico Horizontal Equipment Rail
- Available in 1 or 2 foot lengths

Part Number: ARS-RS042-012

FRS Equipment Roll Stand



- Fairfield compatible Equipment Rail
- Available in 1 or 2 foot lengths

Part Number: FRS-RS042-012

Equipment Roll Stand



- Using a pole clamp, this roll stand can be completely customized for your needs

Part Number: RS1-42

Roll Stand Accessories



Hose Assemblies

- Used to connect a medical gas block to a source
- Custom made to order

Call for part numbers.



Gas Blocks

- Snap-on outlets 1, 2 or 3
- Available for all gases
- All outlet fittings available

Amico and Fairfield gas block options shown



Canisters & Canister Holders

- Disposable and reusable canisters available
- Ring bracket and canister sold separately
- 1200cc, 2000cc, and 3000cc size available

Call for part numbers.



Vacuum Slides

- Conveniently mount your canister suction ring bracket


Part Number: SLIDE2-P (includes two slides)



Blender Bracket

- Used for mounting oxygen/air blenders

Part Number: BLBRKT-R125



Cord Wrap

- Ideal for wrapping wires, hoses and tubing

Part Number: CORDWRAP-R125



Wire Basket

- 10"W x 6"D x 9"H white wire basket
- Convenient storage

Part Number: 100609-R125-BB



Humidifier Bracket


- Used for mounting heated humidifiers

Part Number: HUMBRKT-R125

Please visit www.amico.com for more information on our complete line of mobile solutions.



Distributed By:



SP-RSQ062

the right connection