



Respiratory Equipment



Features:

- Equipment can easily be locked down or moved for security and convenience
- Improves patient safety and care
- Easily mounts to any rail system, boom, or headwall
- Improves workflow and accessibility
- Layout of room can be tailored based on the needs of the provider and patient
- Setup can be customized along the rail to accommodate multiple room configurations
- No assembly required
- 5 year standard warranty

Amico Accessories Instrument Mounts offer unique benefits to healthcare providers including flexibility, adaptability, productivity, aesthetics and ergonomics. Patient rooms can be re-configured quickly and easily in response to changing demographics, protocols or new regulations.




Mounting Options

- A** Amico Horizontal Rail Mount
(Compatible with Amico and Modular)
- V** Amico Vertical Rail Mount
(Compatible with Amico and Modular)
- F** Fairfield Rail Mount (push button optional)
- E** DIN 25mm Rail Mount (30mm, 35mm optional)
- C** 1 ¼" Pole Mount
- M** Monitor Rail Mount
- G** Vertical Rail Mount
(Compatible with Amico, Modular and Hill Rom)
- S** HSI Mount-Vertical Rail Mount

Amico Accessories Instrument Mounts include the adapters required for the above rail systems.

Respiratory Accessories

Resuscitation Bag Holder



- Dual hook design provides an ideal space for most resuscitation bag holders
- Anodized and chrome finishes are optimal for easy cleaning


A RESUSHLD-A

F RESUSHLD-F

V RESUSHLD-V

E RESUSHLD-E25

E-Cylinder Tank Holder



- Holds E-Cylinder tanks
- Available in a white powder coat finish
- Secure locking mechanism


A ECYLHLDR-A

F ECYLHLDR-F

V ECYLHLDR-V

E ECYLHLDR-E25

Blender Bracket



- Anodized Aluminium prevents scratches and slows down wear
- Compatible with all makers of O₂ Blenders


A BLBRKT-A

F BLBRKT-F

V BLBRKT-V

E BLBRKT-E25

Supply Clips



- Ideal for holding single use packaged items
- Quick pull release


A SUPCLIP-2-A

F SUPCLIP-1-F

V SUPCLIP-1-V

E SUPCLIP-1-E25

Humidifier Bracket



- Fisher & Paykel heated humidifier mount
- Bottom support for secure mounting


A HUMBRKT-A

F HUMBRKT-F

V HUMBRKT-V

E HUMBRKT-E25

Tubing Clips



- Neatly organizes O₂ Cannulas and tubes
- Ideal use for Suction Canister tubing
- Useful for hoses, electrical , and nurse call cords


A TUBECLIP-2-A

F TUBECLIP-1-F

V TUBECLIP-1-V

E TUBECLIP-1-E25

Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings and more Respiratory Mounting Solutions.



Distributed By:

