


Wall/Console Slide Bracket & Blank Wall Assembly

General Specifications:

Materials

Slide: Aluminum
Screws: Plated Steel
Cover Plate: Cycloy® (N)

Note: All Wall/Console Slide Bracket and Blank Wall Assemblies are available in M.R.I. versions (see page 2)



Wall Blank (Left) & Wall Slide (Right) shown above



Console Blank (Left) & Console Slide (Right) shown above

Model Numbers

O1-XXXXXXXX-X

The Letter "XXX" Defines the Type:

Wall Outlet = WAL
Console Outlet = CON


The Letter "XXXX" Defines the Type:

Slide = SLIDE
Blank = BLANK

The Letter "X" Defines the Type:

Cycloy = N

Example: Console Blank = O1-CONBLANK-N

Technical Specifications:


M.R.I. Model Numbers

O1-XXXXXXXXM-X

The Letter "XXX" Defines the Type:

Wall Outlet = WAL
Console Outlet = CON

The Letter "XXXXX" Defines the Type:

Slide = SLIDE
Blank = BLANK

M = M.R.I. Compatible

The Letter "X" Defines the Type:

Cycloloy = N

Example: Console Blank (M.R.I.)= O1-CONBLANKM-N